



**Centre of Forensic Sciences
Investigators and Submitters**

**Technical Information Sheet
Toxicology**

October 2024

Contents

Contents	1
Introduction	2
Examination Strategy and Capability.....	2
Urgent Cases	2
Examination.....	2
Instrumentation	3
Glossary.....	4
Appendix 1 – Screening Methods.....	6
GC/MS Screen.....	6
QTOF Screen.....	8
Head-space GC-FID analysis for volatiles (screen and quantitation).....	8
Volatile screen (qualitative only).....	8
Appendix 2 – Drugs Requiring Targeted Analysis.....	9
Appendix 3 – Initial Analyses by Case Type ^a	10
Appendix 4 – No Method Available.....	11
Appendix 5 – Capability of Quantitative Methods	12
Barbiturate method (LC-MS/MS)	12
GHB/BHB method (LC-MS/MS)	12
Cannabinoid method (LC-MS/MS).....	12
LC-MS/MS Mix 2	12
LC-MS/MS Mix 3	12
LC-MS/MS Mix 4	13
LC-MS/MS Mix 5	13
Appendix 6 – Capability of Targeted Qualitative Methods	13

Introduction

The Toxicology Section performs analyses on biological samples (e.g., blood, urine, liver) to determine the absence/presence/concentration(s) of drugs, including alcohol and poisons.

This document is intended as a convenient investigative reference but should not be relied upon as definitive or exhaustive. Please contact the Centre of Forensic Sciences (CFS) Toxicology Section for assistance with questions of an analytical or toxicological nature by e-mail or telephone (647-329-1400 or 647-329-1430). When calling please ask for the appropriate coordinator:

Coroner's Coordinator:

CFSToxicologyCoronerCoordinator@ontario.ca

Criminal Coordinator:

toxcrim@ontario.ca

Examination Strategy and Capability

The screening methods employed in the Toxicology Section are:

1. Gas Chromatography/Mass Spectrometry (GC/MS)
2. Head-Space GC analysis for volatiles
3. Quadrupole Time-of-Flight MS (QTOF)

The targeted/quantitation methods employed in the Toxicology Section are:

1. GC
2. Liquid Chromatography-Mass Spectrometry (LC-MS/MS)
3. Head-Space GC analysis for volatiles

Capabilities of screening methods are presented in Appendix 1. While these screening methods have wide-ranging capabilities, not all drugs may be reliably detected. Appendix 2 contains a list of compounds that may not be identified by the screening methods but may be detected/quantitated by targeted methods. Many of the compounds contained in this list will not be tested for unless specifically requested. If use of a specific drug is known or suspected and is relevant it should be noted in the case synopsis.

The examination strategy, i.e., determining which tests will be performed in a case, is informed by a variety of sources including case type, case history, nature of submitted samples, analytical protocols and capabilities, and discussions with clients. The initial toxicological analyses conducted for a variety of case types are presented in Appendix 3.

Urgent Cases

Requests for expedited analyses must meet specific criteria before being accepted as an urgent case. This process requires authorization by Toxicology Section management.

Examination

All items are visually examined on receipt to check the seal numbers (if present), the contents, and the integrity of the packaging.

Instrumentation

Chromatography: Gas Chromatography (GC); Liquid Chromatography (LC)

Chromatography is an analytical technique used to separate compounds based on their chemical and structural properties. GC uses a pressurized gas, while LC uses a pressurized liquid, in the separation of compounds.

Mass Spectrometry (MS)

MS detects, identifies, and quantitates compounds. An MS can be coupled with a GC or an LC.

Quadrupole Time-of-Flight-MS (QTOF)

QTOF detects and identifies compounds. A QTOF is coupled with an LC.

Tandem MS (MS/MS)

MS/MS detects, identifies, and quantitates compounds and is commonly coupled to a GC or LC.

Ultraviolet and Visible (UV/VIS) Spectrophotometry

UV/VIS spectrophotometry identifies and/or quantitates a drug based on its UV and/or visible light-absorbing properties.

Carbon Monoxide

Carbon monoxide is analyzed by visible spectrophotometry. Results are expressed as % carboxyhemoglobin saturation.

Interpretation

Quantitative results may be expressed as 1) a concentration or 2) as < or > a concentration, e.g., when sufficient for interpretation. Blood ethanol interpretations provided in reports are generally limited to cases in which the detected concentration may be associated with fatalities, may be influenced by post-mortem artefacts, or may have toxic interactions with other drugs.

Measurement Uncertainty

Measurements made with all scientific instruments are associated with variability. No measurement is exact but is an estimate of the true value. Calculation of measurement uncertainty (MU) employs statistical methods to determine the range of values within which the quantitative result is likely to reside. The MU provides a reasonable estimate of the variability associated with the analytical method and is based on the analysis of matrix-matched quality control samples. A minimum of 10 such analyses are used. The MU is calculated with a confidence of 95.45 per cent using a k-factor based on the degrees of freedom as determined by the Student's *t*-test and the standard deviation of the associated quality control data. The MU is expressed in the same units in which the quantitative result is reported, e.g., ng/mL, mg/L and is reported as: quantitative result \pm MU.

Limitations

The focus of this laboratory is drug toxicity. Clinical blood/urine chemistry analysis, e.g., electrolytes, cell counts, gas saturation, creatinine, is not performed. Analysis for antiepileptic drugs is limited to determining drug toxicity, when warranted, based on case history. This laboratory does not have validated methods to analyze some sample types, e.g., oral fluid, hair, bile, muscle, brain tissue. There are a variety of analytical issues that may prevent the detection of some of the drugs that this laboratory is commonly capable of detecting, which include:

- matrix effects
 - degree of putrefaction
 - type of sample (e.g., splenic blood)
 - post-mortem interval
 - storage conditions
- volume of sample submitted
- low concentration of the drug/sensitivity of the method

Conversely, some novel, or rarely encountered, drugs not listed in Appendix 1 may be identified by the GC/MS or QTOF screens. In this case, analytical reference material would be acquired (if available) then analysed to confirm identity. There are drugs/compounds for which the CFS Toxicology Section does not have a method, examples of which are provided in Appendix 4.

Glossary

Abbreviations

Analytical results are reported in terms of mg/100 mL, mg/L, or ng/mL, as shown below:

g	gram
mg	milligram
ng	nanogram
L	litre
mL	millilitre

Breakdown Product

A compound produced either inside or outside the body that may or may not be pharmacologically active.

Carboxyhemoglobin saturation

The percentage of hemoglobin bound by carbon monoxide.

Central Nervous System Depression (CNS depression)

A lowering of the functional activity of the brain and/or spinal cord. Depression of the respiratory and the cardio-regulatory centres are most relevant toxicologically.

Confirmation

The process of verifying the presence of a drug by replicate analysis using the same or different analytical technique(s).

Coroner's Case Analytical Summary

Contains analytical results with the fatal reference and limitations. The Coroner's Case Analytical Summary is accompanied by an Interpretive Guide with information specific to this report type.

Detected

The drug has been identified in the sample. Identification is based on criteria specific to the analytical technique.

Fatal Reference

A minimum drug concentration at which death has been reliably reported in the forensic literature.

Inconclusive

The presence or absence of a drug could not be determined.

Metabolite

The product of enzymatic conversion of a drug within the body to a different compound that may or may not be pharmacologically active.

No [other] significant findings by a [method name(s)]

This comment is inserted to provide a reference to the methods that were used. Appendices 1 and 5 can be used to identify compounds not listed and that were either not detected or the results were deemed to not be toxicologically significant, e.g., caffeine or nicotine. This may also apply to endogenous compounds, e.g., acetone < 2 mg/100 mL.

Not Detected

The drug is either not present or is present but at an amount that cannot be discerned from other constituents in the sample.

Post-mortem redistribution

A phenomenon that refers to a change (either an increase or a decrease) in blood drug concentration after death; post-mortem redistribution may occur regardless of sampling site but is commonly observed as increased drug concentrations in heart blood as compared to femoral blood.

Putrefaction

The decomposition of organic material that involves micro-organisms.

Report

Contains a comprehensive summary of analytical results accompanied by interpretative conclusions.

Therapeutic

The detected drug concentration is generally considered to not be toxicologically significant. The use of this term does not imply clinical efficacy.

Traces

The drug is detected; the concentration is less than the limit of quantitation of a targeted analysis.

Unconfirmed

A drug has been identified by a single procedure but not quantitated or confirmed by a second analysis. Unconfirmed findings may or may not be toxicologically significant.

Appendix 1 – Screening Methods

Drugs that can be reliably detected by screening methods

GC/MS Screen

A	
alpha-pyrrolidinovalerophenone (α-PVP)	cotinine
acetylfentanyl ²	cyclobenzaprine ²
amantadine	cyproheptadine
amitriptyline ²	D
amlodipine ²	desipramine ²
amoxapine ²	dextromethorphan ²
amphetamine ²	dextrorphan*
amphetamine (4-fluoro)	diazepam ²
anabasin	diazepam (nor) ²
atomoxetine	dibucaine
atropine/hyoscyamine	dihydrocodeine
B	diltiazem ²
benzocaine	└ desacetyldiltiazem ²
benzofuran (6-(2-aminopropyl), 6-APB)	dimethyltryptamine
benztropine ²	diphenhydramine ²
benzylpiperazine (BZP) ²	doxepin ²
bromo-dragonfly	doxylamine ²
brompheniramine ²	E
bupivacaine	ephedrine* ²
bupropion ²	estazolam
butylone/ethylone	etizolam ²
butyryl fentanyl ²	ethylone/butylone
C	F
caffeine ²	x-fluoroamphetamine
carbamazepine ¹	fluoxetine ²
cathinone (cath)	fluoxetine (nor) ²
n-ethyl-cath	flurazepam ²
4-fluorometh-cath	flurazepam (n-desalkyl) ²
3-methoxymeth-cath	fluvoxamine ²
4-methyleth-cath	H
meth-cath	haloperidol ²
chlorcyclizine	hydrocodone ²
chlordiazepoxide ²	hydroxychloroquine
chloroquine	hydroxyzine
chlorpheniramine ²	I
chlorpromazine	ibogaine
cisapride	imipramine ²
citalopram* ²	K
clomipramine ²	ketamine ²
clonidine ²	L
clozapine ²	lamotrigine ²
cocaethylene ²	
cocaine ²	
codeine ²	

laudanosiine
levamisole
lidocaine²
loratadine
loxapine²

M

maprotiline
meclizine
mefloquine
meperidine²
meperidine (nor)²
mephedrone²
mepivacaine
methadone²
methamphetamine²
methamphetamine (4-fluoro)
methedrone
methotrimeprazine²
methylenedioxyamphetamine (MDA)²
methylenedioxyethylamphetamine (MDEA)²
methylenedioxymethamphetamine (MDMA)²
3,4-methylenedioxypyrovalerone (MDPV)²
methylone²
methylphenidate²
metoclopramide
metoprolol²
midazolam²
mirtazapine²
moclobemide

N

nicotine²
nortriptyline²

O

olanzapine²
orphenadrine²
oxybutynin
oxycodone²

P

paroxetine²
pentadrone
pentazocine²
pentoxyphylline²
pentylone
phenacetin
phencyclidine (PCP)²
phenethylamines (2C-B, 2C-B-Fly, 2C-T-7, PEA)

pheniramine²
phenmetrazine
phentermine
piperazine, 1-3 chlorophenyl (mCPP)
piperazine, trifluoromethylphenyl (TFMPP)
p-fluorofentanyl
p-methoxyamphetamine (PMA)²
p-methoxymeth-amphetamine (PMMA)
procaine
prochlorperazine²
procyclidine
propoxyphene²
propranolol²
pseudoephedrine*²

Q

quetiapine²
quinidine

R

ropivacaine

S

scopolamine (hyoscine)
sertraline²
strychnine

T

tapentadol
terbinafine
ticlopidine
tramadol²
trazodone²
trihexphenidyl²
trimethoprim
trimipramine²
triprolidine²

V

valeryl fentanyl
varenicline
venlafaxine²
venlafaxine (O-desmethyl)²
verapamil²

X

xylometazoline

Z

zolpidem²
zopiclone breakdown product

QTOF Screen

The QTOF screen is a powerful and sensitive method that can reliably detect the drugs included in the following methods (details are listed in Appendices 5 and 6):

- LC-MS/MS Mix 2
- LC-MS/MS Mix 3 (except carfentanil)
- LC-MS/MS Mix 4
- LC-MS/MS Mix 5 (except: diflunisal, furosemide, ibuprofen, salicylate, vigabatrin)

In addition, the QTOF screen can identify psilocin. The list of drugs potentially identifiable by QTOF is too extensive to list within this document. For questions about a specific drug not listed, please contact the appropriate [case coordinator](#).

*The GC/MS screen and QTOF screen are not capable of distinguishing racemates, therefore compounds such as dextrorphan/levorphanol, citalopram/escitalopram and ephedrine/pseudoephedrine cannot be separated. Similarly, the GC/MS screen cannot distinguish between 2-fluoroamphetamine, 3-fluoroamphetamine, and 4-fluoroamphetamine and the QTOF screen cannot distinguish between methyl fentanyl, butyryl fentanyl, and isobutyryl fentanyl.

Head-space GC-FID analysis for volatiles (screen and quantitation)

acetone
ethanol
isopropanol

methanol
n-propanol (qualitative)

Volatile screen (qualitative only)

difluoroethane
dichloromethane
1,1,1,2-tetrafluoroethane
ethyl acetate
diethyl ether
dimethyl ether

propane
butane
isobutane
toluene
methanol
ethanol

acetone
methyl ethyl ketone
isopropyl alcohol
acetaldehyde
chloroform
gasoline

Appendix 2 – Drugs Requiring Targeted Analysis

Compounds that may not be identified by screening methods but might be detected and/or quantitated by targeted methods.

C

carbon monoxide⁴
cyanide²

D

diflunisal²

F

formic acid³
furosemide²

I

ibuprofen²

T

toluene³

V

valproic acid³
vigabatrin²

Methods used for the quantitation of compounds identified in the preceding appendices are denoted as follows:

¹ GC-NPD

² LC-MS/MS

³ GC-FID

⁴ Visible spectrophotometry

Appendix 3 – Initial Analyses by Case Type^a

Alcohol-impaired driving:	Ethanol
Attempted murder:	dependent upon case history
Confirmation of ketoacidosis:	Ethanol (includes acetone), BHB
Death of child < 5 years of age	Ethanol, QTOF Screen, LC-MS/MS Mix 3, Cannabinoid method
Drug-impaired driving:	QTOF Screen, Cannabinoid method, UDM, GHB
Fatal motor vehicle collision (driver) and aviation death:	Ethanol, QTOF Screen, LC-MS/MS Mix 3, Cannabinoid method, CO ^b
Fire-related death^c:	CO (whole blood required)
Homicide:	Ethanol, QTOF Screen, LC-MS/MS Mix 3, Cannabinoid method
Mandatory inquest:	Ethanol, QTOF Screen, LC-MS/MS Mix 3, Cannabinoid method
Possible drug-related death:	Ethanol, QTOF Screen, LC-MS/MS Mix 3
Rule Out/exclusionary Toxicology:	Ethanol, LC-MS/MS Mix 3
Sexual assault^a:	dependent upon case history
SIU death investigation:	Ethanol, QTOF Screen, LC-MS/MS Mix 3, Cannabinoid method

^a dependent upon sample volume

^b if fire is involved

^c other analyses may be performed dependent upon evidence/suspicion of intoxication

Appendix 4 – No Method Available

Examples of drugs/compounds for which this laboratory does not have a method

Animal toxins

α -bungarotoxin
conotoxin
maurotoxin
tetrodotoxin

Anesthetic gases

halothane
isoflurane
nitrous oxide

Curare-related toxins

alloferine
toxiferine
tubocurarine

Other

insulin
lead, mercury
lithium
polychlorinated biphenyls (PCB)
succinylcholine
thallium

Appendix 5 – Capability of Quantitative Methods

Barbiturate method (LC-MS/MS)

amobarbital (qualitative)
butalbital
pentobarbital
phenobarbital
phenytoin
primidone
secobarbital

Cannabinoid method (LC-MS/MS)

tetrahydrocannabinol (THC)
THC (11-nor-carboxy; THC-COOH)
THC (11-hydroxy; THC-OH, qualitative)
cannabidiol
cannabinol

GHB/BHB method (LC-MS/MS)

γ -hydroxybutyrate (GHB)
 β -hydroxybutyrate (BHB)

LC-MS/MS Mix 2

benztropine
benzylpiperazine
brompheniramine
caffeine (semi-quantitative)
clonidine

ephedrine
haloperidol
ketorolac
loperamide (qualitative)
lidocaine (semi-quantitative)

mitragynine (qualitative)
nicotine (semi-quantitative)
pseudoephedrine
trimeprazine (qualitative)
warfarin

LC-MS/MS Mix 3

6-monoacetylmorphine (6-MAM;
qualitative)
acetyl fentanyl
alprazolam
amitriptyline
amphetamine
benzoylecgonine
bromazolam
bupropion
carfentanil
chlorpheniramine
citalopram/escitalopram
clonazepam
clonazepam (7-amino)
clonazolam
clonazolam (8-amino; qualitative)
cocaethylene
cocaine
codeine
cyclobenzaprine
dextromethorphan
diazepam
diazepam (nor)

diphenhydramine
etizolam
fentanyl
flualprazolam
flubromazolam
flunitrazepam (7-amino)
fluorofentanyl*
fluoxetine
fluoxetine (nor)
flurazepam (n-desalkyl)
hydrocodone
hydromorphone
hydroxyrisperidone/paliperidone
(qualitative)
isotonitazene
ketamine
ketamine (nor)
lysergic acid diethylamide (LSD)
lorazepam
meperidine
meperidine (nor)
mephedrone (qualitative)
methadone

methamphetamine
methylenedioxyamphetamine
methylenedioxyethylamphetamine
methylenedioxymethamphetamine
midazolam
mirtazapine
morphine
nortriptyline
olanzapine
oxazepam
oxycodone
oxymorphone
paroxetine
pseudoephedrine
quetiapine
risperidone
sertraline
temazepam
tramadol (cis)
trazodone
venlafaxine
xylazine
zopiclone

*The method cannot distinguish between para-, meta-, and/or ortho-fluorofentanyl.

LC-MS/MS Mix 4

alprazolam (hydroxyl)	doxylamine	naltrexone
amoxapine	duloxetine	nitrazepam
bromazepam (qualitative)	flunitrazepam	nitrazepam (7-amino)
buprenorphine	flunitrazepam (N-desmethyl)	orphenadrine (qualitative)
butyryl fentanyl	flurazepam	PCP
chlordiazepoxide	fluvoxamine	pentazocine
chlorpromazine	furanyl fentanyl	pheniramine
clobazam	imipramine	promethazine
clomipramine	levorphanol/dextrorphan	propoxyphene
clozapine	(qualitative)	triazolam
demoxepam	loxapine	triazolam (hydroxy)
desipramine	methotrimeprazine	trimipramine
desomorphine	methylenedioxypropylvalerone	U-47700
diltiazem	methylone	venlafaxine (O-desmethyl)
diltiazem (desacetyl)	methylphenidate	ziprasidone
doxepin	naloxone	zolpidem

LC-MS/MS Mix 5

acebutolol	gabapentin	prochlorperazine
acetaminophen	guaifenesin	propafenone
amiodarone	ibuprofen	propranolol
amlodipine	labetalol	pseudoephedrine
atenolol	lamotrigine	salicylate
baclofen	methocarbamol	topiramate
carbamazepine (qualitative)	metoprolol	verapamil
diflunisal	naproxen	vigabatrin
furosemide	pregabalin	

Appendix 6 – Capability of Targeted Qualitative Methods**Urine Drug Mix (UDM; LC-MS/MS)**

6-monoacetylmorphine (6-MAM)	chlordiazepoxide	dextromethorphan
acetyl fentanyl	chlorpheniramine	diazepam
acetyl norfentanyl	citalopram/escitalopram	diazepam (nor)
alprazolam	clobazam	diltiazem
amitriptyline	clomipramine	diltiazem (desacetyl)
amlodipine	clonazepam	diphenhydramine
amoxapine	clonazepam (7-amino)	doxepin
amphetamine	clonazolam	doxylamine
baclofen	clonazolam (8-amino)	duloxetine
benzoylecgonine	clozapine	ephedrine
bromazepam	cocaethylene	etizolam
bromazolam	cocaine	fentanyl
brompheniramine	codeine	fentanyl (nor)
buprenorphine	codeine-6-glucuronide	flualprazolam
buprenorphine glucuronide	cyclobenzaprine	flubromazolam
bupropion	demoxepam	flunitrazepam
butyryl fentanyl	desipramine	flunitrazepam (7-amino)
carfentanil	desomorphine	flunitrazepam (N-desmethyl)

fluoxetine	methamphetamine	pheniramine
fluoxetine (nor)	methylenedioxyamphetamine	pregabalin
flurazepam	methylenedioxyethylamphetamine	propoxyphene
flurazepam (n-desalkyl)	methylenedioxymethamphetamine	propranolol
fluvoxamine	methylenedioxypropylamphetamine	pseudoephedrine
furanyl fentanyl	methylone	quetiapine
gabapentin	methylphenidate	risperidone
GHB	metoprolol	sertraline
heroin	midazolam	tapentadol
hydrocodone	mirtazapine	temazepam
hydromorphone	morphine	temazepam glucuronide
hydromorphone-3-glucuronide	morphine-3-glucuronide	THC-COOH
hydroxyalprazolam	morphine-6-glucuronide	THC-COOH glucuronide
hydroxyrisperidone/paliperidone	naloxone	topiramate
hydroxytriazolam	naltrexone	tramadol (cis)
imipramine	nitrazepam	trazodone
ketamine	nitrazepam (7-amino)	triazolam
ketamine (nor)	nortriptyline	trimipramine
lamotrigine	olanzapine	U-47700
levorphanol/dextrorphan	orphenadrine	venlafaxine
lidocaine	oxazepam	venlafaxine (O-desmethyl)
lorazepam	oxazepam glucuronide	xylazine
lorazepam glucuronide	oxycodone	zaleplon
loxapine	oxymorphone	ziprasidone
meperidine	paroxetine	zolpidem
meperidine (nor)	pentazocine	zopiclone
mephedrone	phenazepam	
methadone	phencyclidine	